

Rollator Nitro SL Blue

Item No. 700300300



- Protrudes with its light weight
- Smartsoft-wheels for a comfortable walking on all undergrounds
- Height adjustable, wide, and padded backrest
- User-friendly height adjustment of handles
- Easy folding mechanism with fold safety
- Can be used indoors and outdoors
- Max. load of bag 5kg
- Available in [blue](#) (700300300), [white](#) (700300200) and [anthracite](#) (700300100)

Technical Data

Order number	700300300
Colour	blue metallic
Overall width	64 cm
Overall length	72 cm
Handle height	83 - 93 cm
Seat width	42 cm
Seat depth	22,5 cm
Seat height	56,5 - cm


Folded width	27 cm
Folded height	83 cm
Folded length	72 cm
Total weight	6.1 kg
Weight without accessories	5.5 kg
Max. Load	136 kg
PU Pal.	20 Pal.
PU pc.	1 Stk.

Rollator Nitro SL Blue

Item No. 700300300



Accessories



Front light, black, for rollators and wheelchairs

Item No.700900014



Lock for rollators and wheelchairs

Item No.700900016



Tray for rollator Nitro / Nitro Carbon / Torro / Torro Acero / Diamond DeLuxe

Item No.724400200



Premium Bag for Nitro SL

Item No.724400137