

# Children's walking aid Nimbo red

Item No. 705100010



- Height-adjustable legs
- Soft rubber wheels provide stability and prevent lateral slippage
- Engageable backstops on the rear wheels prevent rolling backwards
- Front swivel wheels make turning and maneuvering easier
- Front wheels can be locked in one direction
- Latex-free grips with flanged ends
- Easily folded for transport and storage

## Technical Data

|                       |            |
|-----------------------|------------|
| Order number          | 705100010  |
| Colour                | red        |
| Size                  | XS         |
| Body size             | 78 - 98 cm |
| Overall width         | 56 cm      |
| Overall length        | 50 cm      |
| Overall height        | 41 - 51 cm |
| Width between handles | 34,5 cm    |

|                           |         |
|---------------------------|---------|
| Width between rear wheels | 55,5 cm |
| Folded width              | 59 cm   |
| Folded height             | 22 cm   |
| Folded length             | 77 cm   |
| Weight                    | 3.8 kg  |
| Max. Load                 | 30 kg   |
| PU Pal.                   | 30 Pal. |

# Children's walking aid Nimbo red

Item No. 705100010

**drive**

**DeVilbiss**  
HEALTHCARE

## Accessories



**Forearm supports**

Item No.705300200



**Hip / pelvis support**

Item No.705300100



**Seat plate**

Item No.705200100



**Seat harness**

Item No.705300300