

Children's walking aid Nimbo purple

Item No. 705100030



- Height-adjustable legs
- Soft rubber wheels provide stability and prevent lateral slippage
- Engageable backstops on the rear wheels prevent rolling backwards
- Front swivel wheels make turning and maneuvering easier
- Front wheels can be locked in one direction
- Latex-free grips with flanged ends
- Easily folded for transport and storage

Technical Data

Order number	705100030
Colour	purple
Size	XS
Body size	78 - 98 cm
Overall width	56 cm
Overall length	50 cm
Overall height	41 - 51 cm
Width between handles	34,5 cm

Width between rear wheels	55,5 cm
Folded width	59 cm
Folded height	22 cm
Folded length	77 cm
Weight	3.8 kg
Max. Load	30 kg
PU Pal.	30 Pal.

Children's walking aid Nimbo purple

Item No. 705100030



Accessories



Forearm supports

Item No.705300200



Hip / pelvis support

Item No.705300100



Seat plate

Item No.705200100



Seat harness

Item No.705300300