

Children's walking aid Nimbo blue

Item No. 705100120



- Height-adjustable legs
- Soft rubber wheels provide stability and prevent lateral slippage
- Engageable backstops on the rear wheels prevent rolling backwards
- Front swivel wheels make turning and maneuvering easier
- Front wheels can be locked in one direction
- Latex-free grips with flanged ends
- Easily folded for transport and storage

Technical Data

Order number	705100120
Colour	blue
Size	S
Body size	98 - 132 cm
Overall width	58 cm
Overall length	61,5 cm
Overall height	51 - 68 cm
Width between handles	34,5 cm

Width between rear wheels	58 cm
Folded width	62 cm
Folded height	19 cm
Folded length	76 cm
Weight	4.3 kg
Max. Load	50 kg
PU Pal.	30 Pal.

Children's walking aid Nimbo blue

Item No. 705100120

drive

DeVilbiss
HEALTHCARE

Accessories



Forearm supports

Item No.705300200



Hip / pelvis support

Item No.705300100



Seat plate

Item No.705200100



Seat harness

Item No.705300300