

Children's walking aid Nimbo blue

Item No. 705100220



- Height-adjustable legs
- Soft rubber wheels provide stability and prevent lateral slippage
- Engageable backstops on the rear wheels prevent rolling backwards
- Front swivel wheels make turning and maneuvering easier
- Front wheels can be locked in one direction
- Latex-free grips with flanged ends
- Easily folded for transport and storage

Technical Data

Order number	705100220
Colour	blue
Size	M
Body size	116 - 154 cm
Overall width	59 cm
Overall length	73 cm
Overall height	61 - 80,5 cm
Width between handles	38,5 cm

Width between rear wheels	61 cm
Folded width	65 cm
Folded height	17 cm
Folded length	90 cm
Weight	5.4 kg
Max. Load	70 kg
PU Pal.	16 Pal.
PU pc.	1 Stk.

Children's walking aid Nimbo blue

Item No. 705100220

drive

DeVilbiss
HEALTHCARE

Accessories



Forearm supports

Item No.705300200



Hip / pelvis support

Item No.705300100



Seat plate

Item No.705200100



Seat harness

Item No.705300300