

Children's walking aid Nimbo blue

Item No. 705100320



- Height-adjustable legs
- Soft rubber wheels provide stability and prevent lateral slippage
- Engageable backstops on the rear wheels prevent rolling backwards
- Front swivel wheels make turning and maneuvering easier
- Front wheels can be locked in one direction
- Latex-free grips with flanged ends
- Easily folded for transport and storage

Technical Data

Order number	705100320
Colour	blue
Size	L
Body size	150 - 189 cm
Overall width	66,5 cm
Overall length	79,5 cm
Overall height	75,5 - 94,5 cm
Width between handles	41,5 cm

Width between rear wheels	65 cm
Folded width	70 cm
Folded height	20 cm
Folded length	113 cm
Weight	6 kg
Max. Load	90 kg
PU Pal.	14 Pal.
PU pc.	1 Stk.

Children's walking aid Nimbo blue

Item No. 705100320

drive

DeVilbiss
HEALTHCARE

Accessories



Hip / pelvis support

Item No.705300100



Seat plate

Item No.705200100



Seat harness

Item No.705300300



Forearm supports

Item No.705300250