

Children's walking aid Nimbo black

Item No. 705100440



- Height-adjustable legs
- Soft rubber wheels provide stability and prevent lateral slippage
- Engageable backstops on the rear wheels prevent rolling backwards
- Front swivel wheels make turning and maneuvering easier
- Front wheels can be locked in one direction
- Latex-free grips with flanged ends
- Easily folded for transport and storage

Technical Data

Order number	705100440
Colour	black
Size	XL
Body size	184 - 214 cm
Overall width	65 cm
Overall length	93.5 cm
Overall height	94 - 109 cm
Width between handles	48 cm

Width between rear wheels	69,5 cm
Folded width	74 cm
Folded height	21 cm
Folded length	133 cm
Weight	7.8 kg
Max. Load	120 kg
PU Pal.	14 Pal.

Children's walking aid Nimbo black

Item No. 705100440

drive

DeVilbiss
HEALTHCARE

Accessories



Hip / pelvis support

Item No.705300100



Seat harness

Item No.705300350



Seat plate

Item No.705200400



Forearm supports

Item No.705300250