

Positioning Seat Trekker GT blue

Item No. 706000100



- Extensive setting and adjustment option
- Innovative wheel functions such as direction lock, backstop, friction and parking brake are standard on all models
- Stable and light construction ensures safety and versatility
- Tool-free installation of all accessories
- Flexible settings for individual needs and abilities
- Seamless adaptation to developing motor skills
- Motivates kids to move, opens up new spaces and supports therapy work

Technical Data

Order number	706000100
Colour	blue
Size	S
Body size	114 cm
Overall width	57 cm
Overall length	84 cm

Overall height	38 - 47,6 cm
Width between handles	33,5 cm
Width between rear wheels	68 cm
Weight	10.5 kg
Max. Load	34 kg

Positioning Seat Trekker GT blue

Item No. 706000100



Accessories



Forearm support

Item No.TK1035S



Forearm support

Item No.TK1035L



Seat harness

Item No.TK1070S



Seat harness

Item No.TK1070L



Hip belt

Item No.TK1080S



Hip belt

Item No.TK1080M



Hip belt

Item No.TK1080L



Thigh belt

Item No.TK1090S



Thigh belt

Item No.TK1090L



Ankle strap for Trekker

Item No.TK1060S



Ankle strap for Trekker

Item No.TK1060M



Ankle strap for Trekker

Item No.TK1060L



Handle

Item No.TK1045