

Standard rollator Cristallo

Item No. 722100000



- Automatic lock to prevent folding
- Can be used indoors and outdoors
- Multiple reflectors ensure safety
- The folding mechanism in signal color ensures safe and comfortable folding
- The brake can be applied easily and is simple to readjust
- The height of the ergonomic handles is freely adjustable
- Puncture-proof tires
- Max. load of basket and tray 5 kg each
- Supplied as standard with basket, tray, and cane holder

Technical Data

Order number	722100000
Colour	silver
Overall width	61 cm
Overall length	67 cm
Handle height	79 - 100 cm
Seat width	38 cm
Seat depth	16 cm
Seat height	- 61 cm
Width between handles	43 cm
Width between rear wheels	41 cm

Folded width	61 cm
Folded height	86 cm
Folded length	26 cm
Total weight	9.8 kg
Weight	8.8 kg
Weight without basket	9.1 kg
Weight without basket and tray	8,8 kg
Max. Load	130 kg
PU Pal.	16 Pal.

Standard rollator Cristallo

Item No. 722100000



Accessories



Back support for rollators

Item No.710100500



Tilt aid retrofit set

Item No.710100300



Basket

Item No.710000001



Tray

Item No.710000002



Rollator bag

Item No.710000075



Rollator bag

Item No.710000076



One-hand brake

Item No.710000081



Cane holder

Item No.710000012