



- Cushioned and individually adjustable forearm supports
- Automatic lock to prevent folding
- Multiple reflectors ensure safety
- The folding mechanism in signal color ensures safe and comfortable folding
- Dual-function brake (for braking and parking)
- Puncture-proof tires
- Wheel diameter/width 190/50 mm
- Max. load of basket and tray 5 kg each
- Supplied as standard with basket, tray, and cane holder

Technical Data

Order number	722500000
Colour	silver
Overall width	65 cm
Overall length	67 cm
Seat width	38 cm
Seat depth	16 cm
Seat height	- 61 cm
Width between handles	38 cm
Width between rear wheels	43 cm

Höhe Armauflage	102 - 118 cm
Folded width	61 cm
Folded height	70 cm
Folded length	26 cm
Total weight	12.3 kg
Weight	12 kg
Weight without accessories	12 kg
Max. Load	120 kg
PU Pal.	16 Pal.

Accessories



Tilt aid retrofit set

Item No.710100300



Basket

Item No.710000001



Tray

Item No.710000002



Rollator bag

Item No.710000075



Rollator bag

Item No.710000076



Cane holder

Item No.710000012