

Standard rollator Migo 2G

Item No. 723600000



- Extended push handles, now up to 100 cm high
- Tool-free mounting of the cane holder
- 50 mm-wide wheels with closed rim for ease of cleaning
- Can be used indoors and outdoors
- Dual-function brake (for braking and parking)
- Easy to transport – ideal for trips
- Hanging basket and tray easy to install
- Anatomically shaped, non-slip hand grips
- Protection against folding through engaged seat
- Max. load of basket and tray 5 kg each
- Supplied as standard with basket, tray, and cane holder

Technical Data

Order number	723600000
Colour	silver
Overall width	58 cm
Overall length	61 cm
Handle height	80 - 100 cm
Seat width	36 cm
Seat depth	23.5 cm
Seat height	- 60 cm
Width between handles	44,5 cm
Width between rear wheels	43 cm

Folded width	58 cm
Folded height	74 cm
Folded length	23 cm
Total weight	9.1 kg
Weight	8.3 kg
Weight without accessories	8.3 kg
Max. Load	130 kg
PU Pal.	24 Pal.
PU pc.	2 Stk.

Standard rollator Migo 2G

Item No. 723600000



Accessories



Back support for rollators

Item No.710100500



Tilt aid retrofit set

Item No.710100300



Basket

Item No.710000001



Tray

Item No.710000002



Rollator bag

Item No.710000075



Rollator bag

Item No.710000076



One-hand brake

Item No.710000081



Cane holder

Item No.710000012