

Rollator Nitro

Item No. 724100300

drive

DeVilbiss
HEALTHCARE




- **NEW: from now on with soft wheels as standard**
- Can be used indoors and outdoors
- Increased stability and safety thanks to double cross struts
- Extra-large 25 cm steering wheels for optimal steering and pushing comfort
- Passive lighting at the rear, front, and sides as standard
- Cushioned seat surface and backrest as standard
- Simple handle height adjustment
- Cane holder included as standard
- Ergonomic, height-adjustable handles
- Dual-function brake (for braking and parking)
- One-hand folding mechanism ensures safe and comfortable folding
- Max. load of mesh bag 5 kg
- Supplied as standard with cane holder, tilt aid, and modern mesh bag with zipper

Technical Data

Order number	724100300
Colour	champagner
Overall width	52 cm
Overall length	70.5 cm
Handle height	71.8 - 77 cm
Seat width	37 cm
Seat depth	20 cm
Seat height	- 45 cm
Width between handles	34,5 cm

Front wheel diameter	245/33 mm
Rear wheel diameter	200/34 mm
Folded width	25 cm
Folded height	71,8 cm
Folded length	70,5 cm
Total weight	6.7 kg
Max. Load	130 kg
PU Pal.	8 Pal.

Accessories



Tilt aid for Nitro

Item No.724400052



Cup holder for Nitro / Nitro Carbon

Item No.724400085



Tray for rollator Nitro / Nitro Carbon / Torro / Torro Acero / Diamond DeLuxe

Item No.724400200