

Rollator Nitro

Item No. 724200000

drive

DeVilbiss
HEALTHCARE



- **NEW: from now on with soft wheels as standard**
- Can be used indoors and outdoors
- Increased stability and safety thanks to cross struts
- Extra-large 25 cm steering wheels for optimal steering and pushing comfort
- Passive lighting at the rear, front, and sides as standard
- Cushioned seat surface and backrest as standard
- Ergonomic, height-adjustable handles
- Dual-function brake (for braking and parking)
- One-hand folding mechanism ensures safe and comfortable folding
- Max. load of mesh bag 5 kg
- Supplied as standard with cane holder, tilt aid, and modern mesh bag with zipper

Technical Data

Order number	724200000
Colour	champagner
Overall width	61 cm
Overall length	70.5 cm
Handle height	84 - 96.5 cm
Seat width	44.5 cm
Seat depth	22 cm
Seat height	- 57 cm
Width between handles	45 cm

Front wheel diameter	245/33 mm
Rear wheel diameter	200/34 mm
Folded width	25 cm
Folded height	86 cm
Folded length	70,5 cm
Total weight	7.2 kg
Max. Load	130 kg
PU Pal.	8 Pal.

Rollator Nitro

Item No. 724200000



Accessories



Bag

Item No.724400019



Cane holder

Item No.724400040



Cup holder for Nitro / Nitro Carbon

Item No.724400085



**Tray for rollator Nitro / Nitro Carbon /
Torro / Torro Acero / Diamond DeLuxe**

Item No.724400200