

Portable aspiration pump VacuAide 7305

Item No. 7305D-G



- Efficient aspiration made easy
- Thanks to continuous vacuum regulation, the VacuAide can be adapted to any aspiration situation – both at home and in hospital use
- Operation both via the integrated battery (P version) and directly using a power cord (D version)
- 800 ml single-use secretion receptacle included as standard
- Autoclavable 1,200 ml receptacle optionally available
- The secretion receptacles are equipped with an automatic float shutoff function to prevent fluid from overflowing into the device
- Automatic battery charging on connecting to the power grid
- Manufactured in the USA

Technical Data

Order number	7305D-G
Colour	black
Overall width	17.8 cm
Overall length	20,3 cm
Overall height	22.9 cm
Suction power	27 l/min

Air flow max.	27 l/min
Flow rate	27 l/min
Power Supply	50/60 Hz
Vacuum range	80 - 550 mmHg
Total weight	2.9 kg
PU Pal.	33 Pal.