

Portable aspiration pump VacuAide QSU 7314

Item No. 7314P-NE



- Noise-optimized operation with high aspiration power at the same time
- Can be used anywhere as a portable variant (P version) or operated at home (D version)
- Simple operation through on/off switch and vacuum indicator with integrated rotary controller for continuous adjustment of aspiration power
- Available with 800 ml single-use secretion receptacle or autoclavable 1,200 ml receptacle
- Continuous vacuum regulation from 50 to 550 mmHg as well as high aspiration power of 27 l/min
- The carrying bag enables the device to be operated and the aspiration pump incl. accessories to be transported as well as providing a concealable viewing window to check the receptacle content

Technical Data

Order number	7314P-NE
Colour	white-blue
Overall width	20.3 cm
Overall length	21,6 cm
Overall height	21.1 cm
Suction power	27 l/min

Air flow max.	27 l/min
Power Supply	50/60 Hz
Vacuum range	50 - 550 mmHg
Total weight	3 kg
PU Pal.	30 Pal.