

# Portable suction unit VacuAide 7325

Item No. 7325D-U



The VacuAide 7325 secretion suction pump is economical, maintenance free and powerful at the same time. It is ideally suited for the most common suction therapy in the home, in nursing homes as well as special clinical applications. The range includes two versions: the mains-powered 7325D model and the 7325P model, which is equipped with a long-lasting battery.

- Space-saving design, light and powerful
- Tube holder integrated
- Stepless vacuum regulation from 50 to 550 mmHg
- High suction power – up to 27l/min (free flow)
- Long-life lithium-ion battery with 1 year warranty
- Up to 45 minutes running time
- Fast charging time of only 3 hours
- Functional personal carry bag included (P-version)

## Technical Data

|                |            |
|----------------|------------|
| Order number   | 7325D-U    |
| Colour         | white-blue |
| Overall width  | 17.8 cm    |
| Overall length | 18 cm      |
| Overall height | 17.5 cm    |
| Suction power  | 27 l/min   |

|               |               |
|---------------|---------------|
| Air flow max. | 27 l/min      |
| Power Supply  | 60 Hz         |
| Vacuum range  | 50 - 550 mmHg |
| Total weight  | 1.54 kg       |
| PU Pal.       | 33 Pal.       |