

Lightweight rollator Gigo 2G

Item No. 735100200

drive

DeVilbiss
HEALTHCARE



- Back rest now removable and foldable
- Particularly light and attractive aluminum rollator
- One-hand folding mechanism ensures safe and comfortable folding
- Ergonomic handles can be adjusted to 6 height positions
- Max. load of basket and tray 5 kg each
- Supplied as standard with basket, tray, backrest, cane holder and tilt aid

Technical Data

Order number	735100200
Colour	silver
Overall width	61,5 cm
Overall length	61 cm
Handle height	87 - 96 cm
Seat width	34 cm
Seat depth	33.5 cm
Seat height	- 54 cm

Width between handles	40,5 cm
Folded width	61,5 cm
Folded height	100 cm
Folded length	28 cm
Total weight	7.8 kg
Weight without accessories	7 kg
Max. Load	120 kg
PU Pal.	10 Pal.

Lightweight rollator Gigo 2G

Item No. 735100200

drive

DeVilbiss
HEALTHCARE

Accessories



Basket

Item No.710000001



Tray

Item No.710000002



Rollator bag

Item No.710000075



Rollator bag

Item No.710000076



One-hand brake

Item No.710000081



Cane holder Gigo / Gigo 2G

Item No.710000090