

Rollator Torro

Item No. 740500100

drive

DeVilbiss
HEALTHCARE



- Large and soft seat surface provides comfort
- The cushioned back support can be removed
- Ergonomic, height-adjustable handles
- Max. load of the bag 5 kg
- Supplied as standard with bag, back support, and cane holder

Technical Data

Order number	740500100
Colour	red
Overall width	61.5 cm
Overall length	68 cm
Handle height	78 - 94.5 cm
Seat width	42 cm
Seat depth	22.5 cm
Seat height	- 59 cm
Width between handles	44,5 cm

Width between rear wheels	47,5 cm
Folded width	22 cm
Folded height	78 cm
Folded length	68 cm
Total weight	6.9 kg
Weight	6.9 kg
Max. Load	150 kg
PU Pal.	10 Pal.

Accessories



Tilt aid for rollator Torro

Item No.740500114



Mesh bag for rollator Torro & Torro Acero

Item No.740500109



Cane holder for rollator Torro

Item No.740500115



Tray for rollator Nitro / Nitro Carbon / Torro / Torro Acero / Diamond DeLuxe

Item No.724400200