

Rollator Torro Acero

Item No. 740500900

drive

DeVilbiss
HEALTHCARE



- Sturdy steel frame
- Loadable up to 150 kg
- Only 8,1 kg weight
- Large and soft seat surface
- Supplied with net bag and can holder
- Max. Load of bag 5 kg

Technical Data

Order number	740500900
Colour	anthracite
Overall width	64 cm
Overall length	69 cm
Handle height	79,5 - 95 cm
Seat width	42 cm
Seat depth	22,5 cm
Seat height	- 58,5 cm

Width between handles	45,5 cm
Folded width	21 cm
Folded height	79,5 cm
Folded length	69 cm
Weight	8.1 kg
Max. Load	150 kg
PU Pal.	10 Pal.
PU pc.	1 Stk.


Rollator Torro Acero

Item No. 740500900

drive

DeVilbiss
HEALTHCARE

Accessories



Mesh bag for rollator Torro & Torro Acero

Item No.740500109



Bell, black, for rollators and wheelchairs

Item No.700900011



Tray for rollator Nitro / Nitro Carbon / Torro / Torro Acero / Diamond DeLuxe

Item No.724400200



Umbrella Standard

Item No.700900019



Umbrella Reverse

Item No.700900021