

# Rollator Diamond Deluxe

Item No. 745R



- Individual and attractive aluminum rollator
- From rollator to transport wheelchair in just a few steps
- Supplied with swiveling and removable footrests
- Tool-free installation of wheels and footrests
- With folding cloth seat and soft, cushioned backrest
- The removable cloth bag provides extra space for personal items
- Height-adjustable handles
- Easy folding using strap on seat surface
- Space-saving storage and easy transport
- Wheel diameter/width 200/30 mm
- Supplied as standard with cane holder and cloth bag

## Technical Data

|                       |              |
|-----------------------|--------------|
| Order number          | 745R         |
| Colour                | red          |
| Overall width         | 64 cm        |
| Overall length        | 104.2 cm     |
| Handle height         | 77 - 94.5 cm |
| Seat width            | 40 cm        |
| Seat depth            | 19 cm        |
| Seat height           | - 54 cm      |
| Width between handles | 46 cm        |

|                           |         |
|---------------------------|---------|
| Width between rear wheels | 47,5 cm |
| Folded width              | 28 cm   |
| Folded height             | 100 cm  |
| Folded length             | 63 cm   |
| Total weight              | 10.3 kg |
| Weight                    | 10.3 kg |
| Max. Load                 | 136 kg  |
| PU Pal.                   | 18 Pal. |

# Rollator Diamond Deluxe

Item No. 745R



## Accessories




**Rollator bag**

Item No.710000075



**Rollator bag**

Item No.710000076



**Cane holder**

Item No.745002



**Bag**

Item No.745001



**Tray for rollator Nitro / Nitro Carbon / Torro / Torro Acero / Diamond DeLuxe**

Item No.724400200