

Shower chair DS 250

Item No. 771100200



- Comfortable, safe, and stable
- Large seat surface with drainage holes
- Cross strut ensures maximum stability
- Five-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Anodized aluminum tube legs
- Easy to clean

Technical Data

Order number	771100200
Colour	white
Overall width	53.5 - 57 cm
Overall height	73,5 - 83,5 cm
Seat width	43 cm
Seat depth	41 cm
Seat height	41 - 51 cm

Seat height adjustment	5
Total weight	4.5 kg
Weight	4.5 kg
Max. Load	250 kg
PU Pal.	8 Pal.
PU pc.	1 Stk.