

MED AIRE DuoWave

Item No. 820000200

drive

DeVilbiss
HEALTHCARE



Thanks to its continuous manual duo comfort setting on the unit, the MED AIRE DuoWave overlay enables individual adjustment of the cell pressure to the weight of the patient, in both a lying and a sitting position. The integrated pocket in the cover prevents the foam underlay from slipping.

- Individual adjustment of the cell pressure to the weight of the patient, both lying and sitting
- Individually usable thanks to possible static pressure-relief positioning
- Continuous, manual cell pressure adjustment based on a weight scale
- Bi-elastic, quilted, and breathable cover
- Integrated foam pocket in the cover prevents any foam underlay used from slipping
- Fast repair of the overlay is possible thanks to individual air cell exchange
- Can be used with slight to high risk of decubitus and as support in decubitus therapy (cat. 1-4*)

Important

*For people who can change their lying position regularly themselves or with help!

Technical Data

Order number	820000200
Overall width	90 cm
Overall length	200 cm

Overall height	14 cm
User weight	36 - 150 kg
Max. Load	150 kg