



This mattress replacement system combines flexible user settings with maximum comfort. Approved for a user weight of 20 kg to 200 kg, the MED AIRE FlexWave system offers variable alternating cycles and can be used individually thanks to the care, sitting, and static function. The automatic cell pressure adjustment responds to changes in the patient's position; the cell-on-cell design does away with the need for a foam base for long-term lying comfort.

- Sitting & care function plus static pressure-relief positioning
- Alternating cycle variably adjustable from 10 to 25 minutes
- Automatic cell pressure adjustment in case of patient position change
- Manual weight adjustment possible
- With intelligent, automatic start/stop
- Automatic key lock after 30 seconds
- Bi-elastic, quilted, and breathable cover
- Fast repair of the mattress is possible thanks to individual air cell exchange
- Can be used with slight to very high risk of decubitus and as support in decubitus therapy (cat. 1-4*)

Important

*For people who can change their lying position regularly themselves or with help!

Technical Data

Order number	820000300
Overall width	88 cm
Overall length	200 cm
Overall height	20 cm

User weight	20 - 200 kg
Total weight	3.9 kg
Max. Load	200 kg