



- High-quality HR cold foam with special multi-zone, cube-shaped notches ensures ergonomic body adaptation, increased flexibility, and a pleasant bed climate
- With heel comfort zone for specific pressure distribution in the heel area
- Can be used with slight to medium risk of decubitus and as support in decubitus therapy (cat. 1-4*)

Important:

*For people who can change their lying position regularly themselves or with help

Technical Data

Order number	820100110
Overall width	90 cm
Overall height	12 cm
Anzahl der Raumgewicht-Schichten	1

Raumgewicht	40 kg ³
User weight	135 kg
Total weight	10 kg
Max. Load	135 kg