



- Highly-elastic HR cold foam surface for increased lying comfort
- Stabilizing PU foam base with intelligent notching technology ensures improved breathability and a pleasanter bed climate
- With heel comfort zone for specific pressure distribution in the heel area
- Can be used with slight to high risk of decubitus and as support in decubitus therapy (cat. 1-4*)

Important

*For people who can change their lying position regularly themselves or with help!

Technical Data

Order number	820100200
Colour	blue
Overall width	90 cm
Overall height	14 cm
Anzahl der Raumgewicht-Schichten	2

User weight	30 - 140 kg
Total weight	11 kg
Max. Load	140 kg
PU Pal.	10 Pal.