



- Highly-elastic HR cold foam surface for increased lying comfort
- Stabilizing PU foam base with intelligent notching technology ensures improved breathability and a pleasanter bed climate
- With heel comfort zone for specific pressure distribution in the heel area
- Can be used with slight to high risk of decubitus and as support in decubitus therapy (cat. 1-4*)

Important

*For people who can change their lying position regularly themselves or with help!

Technical Data

Order number	820100210
Overall width	90 cm
Overall height	14 cm
Anzahl der Raumgewicht-Schichten	2

User weight	30 - 140 kg
Total weight	11 kg
Max. Load	140 kg
PU Pal.	10 Pal.