



- Highly-elastic HR cold foam surface for increased lying comfort
- Stabilizing PU foam base with intelligent notching technology ensures improved breathability and a pleasanter bed climate
- With heel comfort zone for specific pressure distribution in the heel area
- Can be used with slight to high risk of decubitus and as support in decubitus therapy (cat. 1-4\*)

Important

\*For people who can change their lying position regularly themselves or with help!

## Technical Data

|                                  |           |
|----------------------------------|-----------|
| Order number                     | 820100220 |
| Overall width                    | 100 cm    |
| Overall height                   | 14 cm     |
| Anzahl der Raumgewicht-Schichten | 2         |

|              |             |
|--------------|-------------|
| User weight  | 30 - 140 kg |
| Total weight | 11 kg       |
| Max. Load    | 140 kg      |
| PU Pal.      | 10 Pal.     |