



- Special, very soft foam top layer
- Middle layer consisting of passively acting, viscoelastic special foam
- Stabilizing HR cold foam base
- Also suitable for cachectic patients
- Can be used with slight to very high risk of decubitus and as support in decubitus therapy (cat. 1-4*)
- Care with MED AIRE FlexWave may be sensible in the case of decubitus category 4 in the sacral region

Important

*For people who can change their lying position regularly themselves or with help!

Technical Data

Order number	820100300
Colour	blue
Overall width	90 cm
Overall height	16 cm
Anzahl der Raumgewicht-Schichten	3

Raumgewicht	40 50 50 kg ³
User weight	20 - 135 kg
Total weight	13 kg
Max. Load	135 kg
PU Pal.	10 Pal.