



- Special, very soft foam top layer
- Middle layer consisting of passively acting, viscoelastic special foam
- Stabilizing HR cold foam base
- Also suitable for cachectic patients
- Can be used with slight to very high risk of decubitus and as support in decubitus therapy (cat. 1-4*)
- Care with MED AIRE FlexWave may be sensible in the case of decubitus category 4 in the sacral region

Important

*For people who can change their lying position regularly themselves or with help!

Technical Data

Order number	820100310
Overall width	90 cm
Overall height	16 cm
Anzahl der Raumgewicht-Schichten	3
Raumgewicht	40 50 50 kg ³

User weight	20 - 135 kg
Total weight	13 kg
Max. Load	135 kg
PU Pal.	10 Pal.