



The GRAVIMED Baritexx is a three-layer comfort mattress with good pressure-reducing and pressure-distributing characteristics. It is also suitable for obese patients.

- Lying area made of viscoelastic foam adapts to the body contour and weight; for optimal pressure relief and distribution with low temperature dependency
- Base and middle layer made of stabilizing HR cold foam
- Integrated notching technology for increased lying comfort
- Can be used with slight to very high risk of decubitus and as support in decubitus therapy (cat. 1-4*)

Important

*For people who can change their lying position regularly themselves or with help!

Technical Data

Order number	820100410
Overall width	90 cm
Overall height	16 cm
Anzahl der Raumgewicht-Schichten	3

User weight	80 - 250 kg
Total weight	13 kg
Max. Load	250 kg
PU Pal.	10 Pal.