



The GRAVIMED Baritexx is a three-layer comfort mattress with good pressure-reducing and pressure-distributing characteristics. It is also suitable for obese patients.

- Lying area made of viscoelastic foam adapts to the body contour and weight; for optimal pressure relief and distribution with low temperature dependency
- Base and middle layer made of stabilizing HR cold foam
- Integrated notching technology for increased lying comfort
- Can be used with slight to very high risk of decubitus and as support in decubitus therapy (cat. 1-4*)

Important

*For people who can change their lying position regularly themselves or with help!

Technical Data

| | |
|----------------------------------|-----------|
| Order number | 820100420 |
| Overall width | 100 cm |
| Overall height | 16 cm |
| Anzahl der Raumgewicht-Schichten | 3 |

| | |
|--------------|-------------|
| User weight | 80 - 250 kg |
| Total weight | 13 kg |
| Max. Load | 250 kg |
| PU Pal. | 10 Pal. |