



- Viscoelastic foam top layer for best pressure distribution
- High resilient foam in the middle and base layer, middle layer with intelligent cuts for a good bed climate
- up to a very high decubitus risk and as also supporting a decubitus therapy (category 1-4*)

Important

*for people who can change their lying position themselves or with help regularly!

Technical Data

Order number	820100430
Overall width	120 cm
Overall height	16 cm
Anzahl der Raumgewicht-Schichten	3

User weight	80 - 250 kg
Total weight	13 kg
Max. Load	250 kg
PU Pal.	10 Pal.