

Seat cushion D-Q-SiT Gel

Item No. 820200200

drive

DeVilbiss
HEALTHCARE



- Optimal adaptation to the ischial tuberosity thanks to visco foam with integrated silicone gel
- Reduction of contact pressure through specific pressure distribution
- Can be used with slight to medium risk of decubitus and as support in decubitus therapy (categories 1-4)
- Special D-Q SiT protect cover - for incontinence protection, reduced shear force, and optimal breathability
- Fits on virtually all standard, lightweight, and multi-function wheelchairs as well as care chairs

Technical Data

Order number	820200200
Overall width	40 cm
Overall length	35 cm

Overall height	7 cm
Weight	1 kg
Max. Load	120 kg