

# Seat cushion D-Q-SiT Air 6

Item No. 820200510

**drive**

**DeVilbiss**  
HEALTHCARE



- Individually adjustable, therefore also suitable as a positioning seat cushion
- Can be used with slight to very high risk of decubitus and as support in decubitus therapy (categories 1-4)
- Special D-Q SiT protect cover - for incontinence protection, reduced shear force, and optimal breathability
- Versatile seat cushion, available in 2 heights, optionally with one or two valves
- Fits on virtually all standard, lightweight, and multi-function wheelchairs as well as care chairs

## Technical Data

Order number	820200510
Overall width	40 cm
Overall length	40 cm

Overall height	6 cm
Weight	1.5 kg
Max. Load	150 kg