



The SuppliTo continuous height-adjustable overbed table is intended for use in the home environment and is designed to provide support for homecare. SuppliTo supports the user and patient at the same time and ensures a very high level of patient comfort thanks to various adjustment options. The table top can be adjusted in inclination and offers maximum safety thanks to special fixations.

- Continuous height adjustment
- Virtually tool-free disassembly
- 4 smooth-running rollers, 2 with brakes
- Flexible use in the patient room
- Individual inclination adjustment of the table top with inclination fixation

## Assembly & functions

## Technical Data

Order number	860000001
Colour	white-gray
Overall width	40 cm
Overall length	670 cm
Overall height	71 - 116 cm
tilt adjustment	0 - 90 °

Width tabletop	40 cm
Length tabletop	60 cm
Total weight	8 kg
Weight	8 kg
Max. Load	10 kg