



The SuppliTo twinTop overbed table consists of a tilt-adjustable and a fixed tabletop. A special inclination fixation ensures maximum safety. It supports both the user and the patient, e.g., during homecare activities or food intake in bed. The height of the SuppliTo twinTop can be individually adjusted. Even at higher table settings, the overbed table remains stable.

- Easy to fold up
- Tool-free assembly and disassembly
- 4 smooth-running rollers, 2 with brakes
- Additional table top
- Flexible use in the patient room

## Assembly & functions

## Technical Data

Order number	860000002
Colour	white-gray
Overall width	40 cm
Overall length	82 cm
Overall height	71 - 93 cm
tilt adjustment	0 - 55° in 7 Schritten

Width tabletop	40 cm
Length tabletop	60 und 20 cm
Total weight	9 kg
Weight	9 kg
Max. Load	10 kg