

Standard wheelchair Rotec

Item No. 900100100



- Simple adjustment from short to long armrests
- Continuous adjustment of leg length (39 - 48 cm)
- Footrests adjustable in five steps
- Wheelbase extension possible
- Back cover height-adjustable for optimal sitting position
- Available in eight seat widths
- Handle height adjustable in four steps from 94 to 101.5 cm
- Seat height adjustment from 51 cm to 56 cm possible in three steps
- Powder-coated surfaces
- Chromed hand rims and leg supports
- Removable and swiveling leg supports
- Aluminum footrests with heel straps can be folded up
- 200 x 50 mm PU steering wheels and 24" x 1 3/8" full-rubber drive wheels
- Extra-reinforced front wheel forks
- Reinforced 24" drive wheels
- Quick-release axles as standard for fast wheel installation or removal
- Optionally available with drum brake for escort
- Crash test according to ISO 7176-19: approved for passenger transport in vehicles in combination with all approved passenger and wheelchair restraint systems
- The Rotec with common part principle: gain flexibility and reduce your stock level – with the common part principle: many Rotec replacement parts and accessories are compatible with all Litec 2G, Litec 2G Plus, Freetec, Ecotec 2G, and Multitec wheelchairs

Technical Data

Standard wheelchair Rotec

Item No. 900100100



Order number	900100100
Colour	black-gray
Drum brake	No
Overall width	57 cm
Overall length	106 cm
Overall height	95 - 102 cm
Seat width	38 cm
Seat depth	43 cm
Seat height	51 - 56 cm
Backrest height	42.5 - 50 cm

Width between armrest	0 cm
Footrest adjustment	39 - 48 cm
Angle of footplate	32 °
Drive wheel dimensions	24" x 1 3/8"
Front wheel dimensions	200 x 5 cm
Folded width	28 cm
Total weight	20.2 kg
Weight	20.2 kg
Max. Load	125 kg
PU Pal.	8 Pal.

Accessories



Power knot system Prosafe

Item No.900100781



Height- and length-adjustable armrest

Item No.900300100



Brake lever extension

Item No.900300110



Anti-tilt aid

Item No.900300120




Anti-tilt aid

Item No.900300121



Folding and lockable leg supports

Item No.900300170




Folding and lockable leg supports

Item No.900300171



Folding and lockable leg supports

Item No.900300172



Stump support

Item No.900300180



Stump support

Item No.900300190



Cane holder for wheelchairs

Item No.900300311



Lap belt L

Item No.900300320



Headrest for wheelchairs

Item No.900300330



PU wheel 150 x 35 mm with wheel fork

Item No.900300340

Standard wheelchair Rotec

Item No. 900100100



Oxygen bottle holder for wheelchairs

Item No.900300350



One-hand brake for escort

Item No.900500767



Wheelchair one-hand drive for Rotec and Freetec

Item No.900100778



Wheelchair / scooter ramp

Item No.STDS1097



?Wheelchair tray? universal

Item No.900300307