

Power knot system Prosafe

Item No. 900100781



- Installation with heavy-duty eyelet
- The adapters ensure quick anchoring of a wheelchair in a motor vehicle
- The power knot critically stabilizes the wheelchair design, e.g. in rear-end collisions
- Operator errors when anchoring the wheelchair are impossible
- The lap belt is firmly attached to the adapter and is thus already anatomically correctly worn before setting off
- A diagonal shoulder belt can be worn
- A wheelchair headrest is possible
- Advantages of installation with a heavy-duty eyelet include the increased power reserve upwards and greater flexibility in using various belt systems

Technical Data

| | |
|--------------|-----------|
| Order number | 900100781 |
| Drum brake | No |

| | |
|-----------------------|------|
| Seat height | - cm |
| Width between armrest | 0 cm |