

Standard wheelchair Rotec

Item No. 900150400



- Simple adjustment from short to long armrests
- Continuous adjustment of leg length (39 - 48 cm)
- Footrests adjustable in five steps
- Wheelbase extension possible
- Back cover height-adjustable for optimal sitting position
- Available in eight seat widths
- Handle height adjustable in four steps from 94 to 101.5 cm
- Seat height adjustment from 51 cm to 56 cm possible in three steps
- Powder-coated surfaces
- Chromed hand rims and leg supports
- Removable and swiveling leg supports
- Aluminum footrests with heel straps can be folded up
- 200 x 50 mm PU steering wheels and 24" x 1 3/8" full-rubber drive wheels
- Extra-reinforced front wheel forks
- Reinforced 24" drive wheels
- Quick-release axles as standard for fast wheel installation or removal
- Optionally available with drum brake for escort
- Crash test according to ISO 7176-19: approved for passenger transport in vehicles in combination with all approved passenger and wheelchair restraint systems
- The Rotec with common part principle: gain flexibility and reduce your stock level – with the common part principle: many Rotec replacement parts and accessories are compatible with all Litec 2G, Litec 2G Plus, Freetec, Ecotec 2G, and Multitec wheelchairs

Technical Data

Standard wheelchair Rotec

Item No. 900150400



Order number	900150400
Colour	black-gray
Drum brake	Yes
Overall width	67,5 cm
Overall length	106 cm
Overall height	95 - 102 cm
Seat width	45 cm
Seat depth	43 cm
Seat height	51 - 56 cm
Seat height adjustment	3
Backrest height	42.5 - 50 cm

Armrest height	23 cm
Width between armrest	0 cm
Footrest adjustment	39 - 48 cm
Angle of footplate	32 °
Drive wheel dimensions	24" x 1 3/8"
Front wheel dimensions	200 x 5 cm
Folded width	30,5 cm
Total weight	22.6 kg
Max. Load	125 kg
PU Pal.	8 Pal.

Accessories



Power knot system Prosafe

Item No.900100783




Height- and length-adjustable armrest

Item No.900300100




Brake lever extension

Item No.900300110



Folding and lockable leg supports

Item No.900300170




Folding and lockable leg supports

Item No.900300171



Folding and lockable leg supports

Item No.900300172



Stump support

Item No.900300180



Stump support

Item No.900300190



Cane holder for wheelchairs

Item No.900300311



Lap belt L

Item No.900300320



PU wheel 150 x 35 mm with wheel fork

Item No.900300340



Oxygen bottle holder for wheelchairs

Item No.900300350



Wheelchair one-hand drive for Rotec and Freetec

Item No.900100778