

Standard wheelchair Rotec XL

Item No. 900150900



- Simple and continuous adjustment from short to long armrests
- Three-position seat height adjustment from 51 cm to 56 cm in the standard model for maximum comfort
- Continuous adjustment of leg length (39 - 48 cm)
- Footrests adjustable in five steps
- Wheelbase extension possible
- Back cover height-adjustable for optimal sitting position
- Available in three seat widths: 51, 56, and 61 cm
- Hand grips for the escort are height-adjustable
- With double cross strut for increased safety
- Quick-release axles as standard
- Optionally available with drum brake

Technical Data

Order number	900150900
Colour	black-gray
Drum brake	Yes
Overall width	78 cm
Overall length	107 cm
Overall height	95 - 107.5 cm
Seat width	56 cm
Seat depth	47 cm
Seat height	51 - 56 cm
Backrest height	43.5 - 50 cm

Armrest height	23 cm
Width between armrest	0 cm
Footrest adjustment	39 - 48 cm
Angle of footplate	32 °
Drive wheel dimensions	24" 1 3/8"
Front wheel dimensions	200 x 50 cm
Folded width	30 cm
Total weight	23 kg
Max. Load	190 kg
PU Pal.	8 Pal.

Standard wheelchair Rotec XL

Item No. 900150900



Accessories




Lap belt XL

Item No.900300321



Brake lever extension

Item No.900300110



Stump support

Item No.900300180



Stump support

Item No.900300190



Cane holder for wheelchairs

Item No.900300311



PU wheel 150 x 35 mm with wheel fork

Item No.900300340



Oxygen bottle holder for wheelchairs

Item No.900300350



?Wheelchair tray? universal

Item No.900300309



Height- and length-adjustable armrest

Item No.900300101