

Lightweight wheelchair Freetec

Item No. 900550700



- Simple adjustment from short to long armrests
- Front fork camber adjustment
- Footrests adjustable in five steps from 0 to 32 degrees
- Wheelbase extension possible
- Available in eight seat widths
- Handle height adjustable to four positions
- Back cover depth and height adjustable
- Armrest reinforced with internal steel insert for greater stability
- Heel straps provide additional hold
- 200 x 50 mm PU steering wheels and 24" x 1 3/8" full-rubber drive wheels
- Extra-reinforced front wheel forks
- Three-position seat height adjustment from 45 cm to 50 cm in the standard model for maximum comfort
- Continuous adjustment of leg length
- Quick-release axles as standard for fast wheel installation or removal
- Optionally available with drum brake for escort
- Crash test according to ISO 7176-19: approved for passenger transport in vehicles in combination with all approved passenger and wheelchair restraint systems
- The Freetec with common part principle: gain flexibility and reduce your stock level – with the common part principle: many Freetec replacement parts and accessories are compatible with all Litec 2G, Litec 2G Plus, Rotec, Ecotec 2G, and Multitec wheelchairs

Technical Data

Lightweight wheelchair Freetec

Item No. 900550700



Order number	900550700	Backrest height	38.5 - 46 cm
Colour	black-gray	Armrest height	21,5 cm
Drum brake	Yes	Width between armrest	0 cm
Overall width	72,5 cm	Footrest adjustment	34.5 - 48.5 cm
Overall length	103 cm	Drive wheel dimensions	24" x 1 3/8"
Overall height	92 - 99.5 cm	Front wheel dimensions	200 x 5 cm
Seat width	52 cm	Folded width	30,5 cm
Seat depth	43 - 50 cm	Total weight	19 kg
Seat height	45 - 50 cm	Max. Load	125 kg
Seat height adjustment	4	PU Pal.	8 Pal.

Lightweight wheelchair Freetec

Item No. 900550700



Accessories



Power knot system Prosafe

Item No.900100783




Height- and length-adjustable armrest

Item No.900300100




Brake lever extension

Item No.900300110



Folding and lockable leg supports

Item No.900300170




Folding and lockable leg supports

Item No.900300171



Folding and lockable leg supports

Item No.900300172



Stump support

Item No.900300180



Stump support

Item No.900300190




Wheelchair tray for Rotec/Freetec

Item No.900300306



Cane holder for wheelchairs

Item No.900300311



Lap belt L

Item No.900300320



PU wheel 150 x 35 mm with wheel fork

Item No.900300340



Oxygen bottle holder for wheelchairs

Item No.900300350



**Wheelbase extension adapter block
Freetec**

Item No.900300360

Lightweight wheelchair Freetec

Item No. 900550700

drive

DeVilbiss
HEALTHCARE



**Wheelchair one-hand drive for Rotec
and Freetec**

Item No.900100778



Seat cushion

Item No.900300272