

# Transport wheelchair D-Lite 12.5

Item No. 900800100

**drive**

**DeVilbiss**  
HEALTHCARE



- Lightweight aluminum frame
- Fold-down backrest
- 12.5" wheels
- Height-adjustable armrests
- Removable and swiveling leg supports
- Available in two seat widths
- Including brake for escort

## Technical Data

|                |              |
|----------------|--------------|
| Order number   | 900800100    |
| Colour         | black-red    |
| Overall width  | 61 cm        |
| Overall length | 89 cm        |
| Overall height | 92.5 - 98 cm |
| Seat width     | 45 cm        |
| Seat depth     | 39.5 cm      |
| Seat height    | - 49 cm      |

|                 |         |
|-----------------|---------|
| Steering wheels | 7"      |
| Folded height   | 73 cm   |
| Folded width    | 28 cm   |
| Folded length   | 67,5 cm |
| Total weight    | 13.2 kg |
| Max. Load       | 125 kg  |
| PU Pal.         | 8 Pal.  |