

Transport wheelchair D-Lite 24

Item No. 900800300



- Lightweight aluminum frame
- Fold-down backrest
- 24" wheels
- Height-adjustable armrests
- Removable and swiveling leg supports
- Available in two seat widths

Technical Data

Order number	900800300
Colour	black-red
Overall width	67 cm
Overall length	101 cm
Overall height	93 - 98.5 cm
Seat width	45 cm
Seat depth	39.5 cm

Seat height	- 49 cm
Folded height	73 cm
Folded width	32 cm
Folded length	69 cm
Total weight	14.4 kg
Max. Load	125 kg
PU Pal.	6 Pal.