

Multifunction wheelchair Multitec

Item No. 900950200



- Adjustable-angle seat unit continuously supported by gas pressure damper
- Adjustable-angle, ergonomically shaped backrest continuously supported by gas pressure damper independently of the seat unit
- Continuous seat depth adjustment
- Specially shaped surface structure, optimal pressure distribution, and stable positioning of the user
- Height-adjustable push handles
- Folding leg supports with knee and calf cushions (height and depth adjustable)
- Adjustable-angle footrests with heel straps
- Anti-tilt supports (height-adjustable) and headrest (optimally adjustable height and inclination) fitted as standard
- Puncture-proof 24" PU drive wheels with quick-release axles and chrome-plated hand rims
- Optionally available with drum brake for escort
- Accessories: abductor wedge, side supports, wheelchair tray, adapter plates to widen seat, and adapter set for pushing aid PowerStroll
- Care aid number: 50.45.07.2033
- Crash test according to ISO 7176-19: approved for passenger transport in vehicles in combination with all approved passenger and wheelchair restraint systems
- The Multitec with common part principle: gain flexibility and reduce your stock level – with the common part principle: many Multitec replacement parts and accessories are compatible with all Litec 2G, Litec 2G Plus, Rotec, Ecotec 2G, and Freetec wheelchairs

Technical Data

Multifunction wheelchair Multitec

Item No. 900950200



Order number	900950200	Backrest angle	-3 - 36 °
Colour	black-gray	Height adjustment armrest	3
Drum brake	No	Armrest height	17 - 20.5 cm
Overall width	65 cm	Width between armrest	0 cm
Overall length	112 - 152 cm	Footrest adjustment	40 - 60 cm
Overall height	97 - 114 cm	Drive wheel dimensions	24"
Seat width	44 cm	Front wheel dimensions	190 cm
Seat depth	42 - 48 cm	Total weight	32 kg
Seat height	- 45 cm	Weight	32 kg
	-1 - 19 °	Max. Load	135 kg
Backrest height	51 - 65 cm	PU Pal.	3 Pal.

Multifunction wheelchair Multitec

Item No. 900950200

drive

DeVilbiss
HEALTHCARE

Accessories



Lap belt XL

Item No.900300321



Rollator bag

Item No.710000075



Rollator bag

Item No.710000076



Folding and lockable leg supports

Item No.900300170



Folding and lockable leg supports

Item No.900300171



Folding and lockable leg supports

Item No.900300172




Stump support

Item No.900300180



Stump support

Item No.900300190



Wheelchair / scooter ramp

Item No.STDS1097



Thermal wheelchair cushion

Item No.RT-CU009



Handmuff

Item No.RT-Handmuff



Abductor wedge

Item No.900600902



Side supports (pair)

Item No.900600900



Wheelchair tray

Item No.900600901

Multifunction wheelchair Multitec

Item No. 900950200

drive

DeVilbiss
HEALTHCARE



Adapter set for PowerStroll

Item No.900600910




Adapter plate for armrest widening

Item No.900900901



Multitec accessory package

Item No.900600903



Seat cushion

Item No.900300271