

Standard wheelchair Ecotec 2G

Item No. 910350300



- Very light and versatile wheelchair
- Armrests can be swiveled away
- Folding footrests with heel straps
- Upholstered seat back and seat cushion
- Height of aluminum wheel forks can be adjusted to three positions
- Powder-coated steel hand rims
- Seat height adjustment of 46 and 51 cm
- Quick-release axles as standard for fast wheel installation or removal
- 175 x 45 mm PU steering wheels and 24" full rubber drive wheels
- Integrated brake lever extension can be extended and folded over
- Crash test according to ISO 7176-19: approved for passenger transport in vehicles in combination with all approved passenger and wheelchair restraint systems
- The Ecotec 2G with common part principle: gain flexibility and reduce your stock level – with the common part principle: many Ecotec 2G replacement parts and accessories are compatible with all Litec 2G, Litec 2G Plus, Freetec, Rotec, and Multitec wheelchairs

Technical Data

Order number	910350300
Drum brake	Yes
Overall width	67 cm
Overall length	107 cm
Overall height	90 - 95 cm
Seat width	46 cm
Seat depth	42 cm
Seat height	49 / 54 cm
Seat height adjustment	2

2 °

Backrest height	39 cm
Armrest height	22.5 cm
Footrest adjustment	36 - 47 cm
Drive wheel dimensions	24"
Front wheel dimensions	7 x 2" cm
Folded width	32 cm
Total weight	20 kg
Max. Load	130 kg
PU Pal.	8 Pal.

Standard wheelchair Ecotec 2G

Item No. 910350300




Accessories



Lap belt XL

Item No.900300321



Infusion holder for wheelchairs

Item No.STDS820



Rollator bag

Item No.710000075



Rollator bag

Item No.710000076



Height- and length-adjustable armrest

Item No.900300100




Anti-tilt aid

Item No.900300120



Anti-tilt aid

Item No.900300121




Folding and lockable leg supports

Item No.900300170



Folding and lockable leg supports

Item No.900300171



Folding and lockable leg supports

Item No.900300172



Stump support

Item No.900300180



Stump support

Item No.900300190



Cane holder for wheelchairs

Item No.900300311



Headrest for wheelchairs

Item No.900300330

Standard wheelchair Ecotec 2G

Item No. 910350300



PU wheel 150 x 35 mm with wheel fork

Item No.900300340




Oxygen bottle holder for wheelchairs

Item No.900300350



One-hand brake for escort

Item No.900500767



Wheelchair / scooter ramp

Item No.STDS1097



One-hand parking brake for wheelchairs

Item No.900500750



One-hand parking brake for wheelchairs

Item No.900500751



Thermal wheelchair cushion

Item No.RT-CU009



Handmuff

Item No.RT-Handmuff



Seat cushion

Item No.900300270