

Lightweight wheelchair Litec 2G

Item No. 910500300



- Very light and versatile wheelchair
- Armrests can be swiveled away
- Adjustable-angle (0° - 30°) and folding footrests as standard
- Leg support height adjustable to five positions
- Upholstered seat back and seat cushion
- Reinforced aluminum frame for optimal stability
- Height of aluminum wheel forks can be adjusted to three positions
- Three-position seat height adjustment from 46 to 52 cm
- Wheelbase extension possible
- 175 x 45 mm PU steering wheels and 24" full rubber drive wheels
- Quick-release axles as standard for fast wheel installation or removal
- Crash test according to ISO 7176-19: approved for passenger transport in vehicles in combination with all approved passenger and wheelchair restraint systems
- The Litec 2G with common part principle: gain flexibility and reduce your stock level – with the common part principle: many Litec 2G replacement parts and accessories are compatible with all Ecotec 2G, Litec 2G Plus, Freetec, Rotec, and Multitec wheelchairs

Technical Data

Lightweight wheelchair Litec 2G

Item No. 910500300



Order number	910500300
Colour	black-gray
Overall width without drum brake	64 cm
Drum brake	No
Overall width	64 cm
Overall length	104 cm
Overall height	87 - 92 cm
Seat width	46 cm
Seat depth	42 cm
Seat height	47 - 52 cm
Backrest height	42 cm

Backrest angle	4 °
Armrest height	23.5 cm
Width between armrest	0 cm
Footrest adjustment	36 - 44 cm
Drive wheel dimensions	24" x 1 3/8"
Front wheel dimensions	190 x 38 cm
Folded width	30 cm
Total weight	14.4 kg
Max. Load	125 kg
PU Pal.	8 Pal.

Lightweight wheelchair Litec 2G

Item No. 910500300

drive


DeVilbiss
HEALTHCARE

Accessories



Lap belt XL

Item No.900300321



Infusion holder for wheelchairs

Item No.STDS820



Rollator bag

Item No.710000075



Rollator bag

Item No.710000076



Height- and length-adjustable armrest

Item No.900300100



Anti-tilt aid

Item No.900300120



Anti-tilt aid

Item No.900300121



Stump support

Item No.900300180



Stump support

Item No.900300190



Cane holder for wheelchairs

Item No.900300311



Headrest for wheelchairs

Item No.900300330



Oxygen bottle holder for wheelchairs

Item No.900300350



One-hand brake for escort

Item No.900500767



Wheelchair / scooter ramp

Item No.STDS1097

Lightweight wheelchair Litec 2G

Item No. 910500300



Thermal wheelchair cushion

Item No.RT-CU009



Handmuff

Item No.RT-Handmuff



Seat cushion

Item No.900300270



Seat cushion

Item No.900300271



Seat cushion

Item No.900300272