

Lightweight wheelchair Litec 2G Plus

Item No. 910600200



- Very light and versatile wheelchair
- Armrests can be swiveled away
- Adjustable-angle (0° - 30°) and folding footrests as standard
- Leg support height adjustable to five positions
- Upholstered seat back and seat cushion
- Reinforced aluminum frame for optimal stability
- Height of aluminum wheel forks can be adjusted to three positions
- Three-position seat height adjustment from 45 to 51 cm
- Wheelbase extension possible
- 175 x 45 mm PU steering wheels and 24" full rubber drive wheels
- Quick-release axles as standard for fast wheel installation or removal
- With premium equipment, including: seat cushion, height- and length-adjustable armrests, cane holder, lap belt & anti-tilt support
- Extensive adjustment options thanks to flexible, adaptable backrest, variable seat depths and push handles
- The Litec 2G Plus with common part principle: gain flexibility and reduce your stock level – with the common part principle: many Litec 2G Plus replacement parts and accessories are compatible with all Litec 2G, Ecotec 2G, Freetec, Rotec, and Multitec wheelchairs

Technical Data

Lightweight wheelchair Litec 2G Plus

Item No. 910600200



Order number	910600200	Backrest angle	4 °
Overall width without drum brake	60 cm	Height adjustment armrest	8
Drum brake	No	Armrest height	25.5 - 32.5 cm
Overall width	60 cm	Footrest adjustment	43 - 51 cm
Overall length	106.5 cm	Drive wheel dimensions	24" x 1 3/8"
Overall height	89.5 - 98.5 cm	Front wheel dimensions	200 x 38 cm
Seat width	42 cm	Folded width	30 cm
Seat depth	42 - 47 cm	Total weight	16.2 kg
Seat height	46 - 51 cm	Max. Load	125 kg
Backrest height	42 - 46 cm	PU Pal.	8 Pal.


Lightweight wheelchair Litec 2G Plus

Item No. 910600200

drive

DeVilbiss
HEALTHCARE

Accessories



Infusion holder for wheelchairs

Item No.STDS820




Rollator bag

Item No.710000075




Rollator bag

Item No.710000076




Folding and lockable leg supports

Item No.900300170




Folding and lockable leg supports

Item No.900300171




Folding and lockable leg supports

Item No.900300172



Stump support

Item No.900300180



Stump support

Item No.900300190



Headrest for wheelchairs

Item No.900300330



Wheelchair / scooter ramp

Item No.STDS1097



Thermal wheelchair cushion

Item No.RT-CU009



Handmuff

Item No.RT-Handmuff



Seat cushion

Item No.900300270