

# Lightweight wheelchair Litec 2G Plus

Item No. 910600400



- Very light and versatile wheelchair
- Armrests can be swiveled away
- Adjustable-angle (0° - 30°) and folding footrests as standard
- Leg support height adjustable to five positions
- Upholstered seat back and seat cushion
- Reinforced aluminum frame for optimal stability
- Height of aluminum wheel forks can be adjusted to three positions
- Three-position seat height adjustment from 45 to 51 cm
- Wheelbase extension possible
- 175 x 45 mm PU steering wheels and 24" full rubber drive wheels
- Quick-release axles as standard for fast wheel installation or removal
- With premium equipment, including: seat cushion, height- and length-adjustable armrests, cane holder, lap belt & anti-tilt support
- Extensive adjustment options thanks to flexible, adaptable backrest, variable seat depths and push handles
- The Litec 2G Plus with common part principle: gain flexibility and reduce your stock level – with the common part principle: many Litec 2G Plus replacement parts and accessories are compatible with all Litec 2G, Ecotec 2G, Freetec, Rotec, and Multitec wheelchairs

## Technical Data

# Lightweight wheelchair Litec 2G Plus

Item No. 910600400



|                 |                |
|-----------------|----------------|
| Order number    | 910600400      |
| Drum brake      | No             |
| Overall width   | 68 cm          |
| Overall length  | 106.5 cm       |
| Overall height  | 89.5 - 98.5 cm |
| Seat width      | 50 cm          |
| Seat depth      | 42 - 47 cm     |
| Seat height     | 46 - 51 cm     |
| Backrest height | 42 - 46 cm     |
| Backrest angle  | 4 °            |

|                           |                |
|---------------------------|----------------|
| Height adjustment armrest | 8              |
| Armrest height            | 25.5 - 32.5 cm |
| Footrest adjustment       | 43 - 51 cm     |
| Drive wheel dimensions    | 24" x 1 3/8"   |
| Front wheel dimensions    | 200 x 38 cm    |
| Folded width              | 30 cm          |
| Total weight              | 16.4 kg        |
| Max. Load                 | 125 kg         |
| PU Pal.                   | 8 Pal.         |

# Lightweight wheelchair Litec 2G Plus

Item No. 910600400

**drive**


**DeVilbiss**  
HEALTHCARE

## Accessories



**Lap belt XL**

Item No.900300321



**Infusion holder for wheelchairs**

Item No.STDS820



**Rollator bag**

Item No.710000075




**Rollator bag**

Item No.710000076



**Folding and lockable leg supports**

Item No.900300170



**Folding and lockable leg supports**

Item No.900300171



**Folding and lockable leg supports**

Item No.900300172



**Stump support**

Item No.900300180



**Stump support**

Item No.900300190



**Headrest for wheelchairs**

Item No.900300330



**Oxygen bottle holder for wheelchairs**

Item No.900300350



**One-hand brake for escort**

Item No.900500767



**Wheelchair / scooter ramp**

Item No.STDS1097



**Thermal wheelchair cushion**

Item No.RT-CU009

# Lightweight wheelchair Litec 2G Plus

Item No. 910600400

**drive**

**DeVilbiss**  
HEALTHCARE



**Handmuff**

Item No. RT-Handmuff



**Seat cushion**

Item No. 900300272