

Lightweight wheelchair Litec 2G Plus

Item No. 910650400



- Very light and versatile wheelchair
- Armrests can be swiveled away
- Adjustable-angle (0° - 30°) and folding footrests as standard
- Leg support height adjustable to five positions
- Upholstered seat back and seat cushion
- Reinforced aluminum frame for optimal stability
- Height of aluminum wheel forks can be adjusted to three positions
- Three-position seat height adjustment from 45 to 51 cm
- Wheelbase extension possible
- 175 x 45 mm PU steering wheels and 24" full rubber drive wheels
- Quick-release axles as standard for fast wheel installation or removal
- With premium equipment, including: seat cushion, height- and length-adjustable armrests, cane holder, lap belt & anti-tilt support
- Extensive adjustment options thanks to flexible, adaptable backrest, variable seat depths and push handles
- The Litec 2G Plus with common part principle: gain flexibility and reduce your stock level – with the common part principle: many Litec 2G Plus replacement parts and accessories are compatible with all Litec 2G, Ecotec 2G, Freetec, Rotec, and Multitec wheelchairs

Technical Data

Order number	910650400
Drum brake	Yes
Overall width	70 cm
Overall length	106.5 cm
Overall height	89.5 - 98.5 cm
Seat width	50 cm
Seat depth	42 - 47 cm
Seat height	46 - 51 cm
Backrest height	42 - 46 cm

Backrest angle	4 °
Height adjustment armrest	8
Armrest height	25.5 - 32.5 cm
Footrest adjustment	43 - 51 cm
Drive wheel dimensions	24" x 1 3/8"
Front wheel dimensions	200 x 38 cm
Folded width	32 cm
Total weight	17.7 kg
Max. Load	125 kg

Lightweight wheelchair Litec 2G Plus

Item No. 910650400

drive

DeVilbiss
HEALTHCARE

Accessories



Infusion holder for wheelchairs

Item No.STDS820




Rollator bag

Item No.710000075




Rollator bag

Item No.710000076



Folding and lockable leg supports

Item No.900300170




Folding and lockable leg supports

Item No.900300171



Folding and lockable leg supports

Item No.900300172



Stump support

Item No.900300180



Stump support

Item No.900300190



Headrest for wheelchairs

Item No.900300330



Oxygen bottle holder for wheelchairs

Item No.900300350



One-hand brake for escort

Item No.900500767



Wheelchair / scooter ramp

Item No.STDS1097



Thermal wheelchair cushion

Item No.RT-CU009



Handmuff

Item No.RT-Handmuff

Lightweight wheelchair Litec 2G Plus

Item No. 910650400

drive

DeVilbiss
HEALTHCARE



Seat cushion

Item No.900300272