

TOURER PLUS (ST4D 2G Plus)

Item No. 950100326



**Please note the changed item codes of our scooter portfolio for the international market. These change as follows:
Item No. 950100336-MD**

If you have any further questions, please contact our customer service.

This scooter is impressive not only because of the sophisticated details but also because of its powerful motor performance. Its many features ensure that you can flexibly arrange your day with the Tourer Plus.

- The latest lighting technology reduces battery consumption
- Approved for transportation in designated buses in the public transport system with pictogram attached
- With the high-low mode, the speed can quickly be reduced to half the maximum speed by pressing a single button
- Identical model also available as 12 km/h scooter
- Lockable steering column which can be adjusted to suit
- Sleek, agile design

Technical Data

Order number	950100326
Colour	red
Overall width	64 cm
Overall length	120 cm
Overall height	115.5 cm
Seat width	46.5 cm
Seat depth	44 - 49 cm
Seat height	- cm
Width between armrest	43 - 55 cm
Speed	12 km/h

Battery capacity 2x	50 Ah
Motor capacity	700 W
Continuous-power	700 W
Gradeability	10° / 17 %
Clearance	6,5 cm
Front tires	290 mm
Rear tires	290 mm
Weight without battery	71,5 kg
Max. Load	136 kg

TOURER PLUS (ST4D 2G Plus)

Item No. 950100326

drive

DeVilbiss
HEALTHCARE

Accessories



Wheelchair / scooter ramp

Item No.STDS1097



Cane holder

Item No.950000028



Rain poncho

Item No.RT-SC001



Rain poncho

Item No.RT-SC004



Scooter garage, folding

Item No.SCGARAGE




Scooter bag for the backrest

Item No.RT-5026



Scooter saddle bag

Item No.RT-5022



Scooter accessory package

Item No.RT-SCPACK



Walking aid/rollator holder

Item No.950000089



Cane holder

Item No.950000029



Scooter garage

Item No.RT-4444